

[Scanningelectronmicroscope for 3dimensional software]

JSM-7000Series MPSeries




Simple and youthful.

- simple
- Easy Operation
- Operation simpleGUI
- \Diamond Quick response
- youthful
- · A variety of data output function
 - Customization feature


(cobservation screen in three steps operation)

(3D display screen)





《observation observation screen in 3 steps operation》《Cross-sectional observation》

image loading 33DConstruction+3Dview ②Alignment



Possible cross-sectional height analysis of any direction







It can be processed in the Excel data measurement result

《A variety of three-dimensional display》


Gray scale display




Color map display



Mesh display




Stereoscopic



《3D measurement model》

Height data created by the 3D model display can measure various



《Operating environment》

- OS: WindowsXP, 7,8 (Japanese / English)
- display: XGA (1024×768) 0ver
- Color Settings : 32-bit color
- Other: recommended wheel mouse
- Supported image formats : BMP / JPEG / TIFF
- corresponding image resolution : maximum of 2560 x 1920

System in Frontier Inc.

 ⊤ 190-0012 Tokyo, 2-8-3
 Shinsuzuharu Building 4F
 Phone number 042-526-4364
 FAX number 042-526-4370
 URL http://www.sifi.co.jp

Tachikawa

Akebono-cho

ik

